

# PRENATAL HEARTBEAT

Fetal Doppler  
Model: 210

The Nidek Fetal Doppler is a hand-held ultrasound transducer used to detect the heartbeat of a fetus for prenatal care.

A Fetal Doppler, sometimes also called a baby Doppler, is used to detect Fetal Heart Rate (FHR) after about 12 weeks pregnancy. It can be used in hospitals, clinics or at home.



## Features

- Pocket Doppler, waterproof ultrasonic probe
- LCD display of heart rate
- Backlight, Built-in speaker
- Audio output can be connected to recorder
- Automatic switch-off
- Three FHR modes: real-time, average, manual
- Battery energy Indicator
- Sound alarm function when FHR is not in the normal range
- International standard, robust construction, ease-of-operation, low maintenance

## Specification

- Ultrasonic Frequency: 2.5 MHz
- LCD size: 37 x 23mm
- FHR display range: 50bpm - 210bpm
- Dimensions: 97mm (Wide) x 37mm (Deep) x 132mm (Height)
- Weight: 0.35kg
- Continuous working time:  $\geq$  8 hours
- Standard accessories: 3.6 V Rechargeable Battery, Gel & Carry bag



### Nidek Medical India Private Limited

**Head Office:** 238B A.J.C. Bose Road, Kolkata 700020  
Tel: (033) 4017 7000 (30 Lines), Fax: 91 (033) 2289 4217

**Delhi:** 238 Ansal Chambers II, 6, Bhikaji Cama Palace, New Delhi 110066  
Tel: (011) 2616 5742, 93501 83222, 93135 77444

**Chennai:** AE - 32, 7th Street, 10th Main Road, Anna Nagar, Chennai 600040  
Tel: (044) 4266 3662, 93808 06777, 93459 70044

**Mumbai:** 211 B, Damji-Shamji Industrial Estate, 9 LBS Marg, Kurla (W), Mumbai 400070  
Tel: (022) 4215 5441 / 93204 22266

**Hyderabad:** H.N0-3-41-66, Sri Laxmi Nilayam, South Block, West Marredpally Secunderabad 500026  
Tel: 93468 32266

**Nagpur:** 93738 38883, **Ahmedabad:** 93284 22333, **Lucknow:** 93361 20941

E-mail: [sales@nidekindia.com](mailto:sales@nidekindia.com), Website: [www.nidekindia.com](http://www.nidekindia.com)