

Sleep Uninterrupted

Auto-CPAP

The NiDEK Auto-CPAP delivers great value, simplicity and quiet performance in a very small package. The NiDEK Auto-CPAP offers an Auto mode with automatic altitude adjustment, in addition to compliance and efficiency monitoring with facility to download data to a mini SD card.


The NiDEK Auto-CPAP flow sensing technology utilizes a special honeycomb structure to create a constant flow resistance for minimizing pressure swings. The net result is a quiet and smooth performance that provides effective Auto-CPAP therapy.

The NiDEK Auto-CPAP also offers options such as auto start, alarm and A-RELIEF technology.



Mode	CPAP, ACPAP
Max Pressure	20 cm H ₂ O
Min Pressure	4 cm H ₂ O
Ramp	0-60 min
Alarm	Power Failure, Mask, Snoring
A-Relief	1-3
Display	2.8 inch Colour
Noise Level	<30db@10 cm H ₂ O
Memory	Micro SD
Dimension	21 x 11 x 11 cm
Humidifier	Optional

Warranty: Two Years



Nidek Medical India Private Limited

Head Office: 238B A.J.C. Bose Road, Kolkata 700020,
Tel: (033) 4017 7000 (30 Lines), Fax: 91 (033) 2289 4217

Delhi: 238 Ansal Chambers II, 6, Bhikaji Cama Palace
New Delhi 110066, Tel: (011) 2616 5742, 93501 83222, 93135 77444

Chennai: AB-17, 5th Street, Anna Nagar, Chennai 600040
Tel: (044) 4266 3662, 93808 06777, 93459 70044

Mumbai: 211 B, Damji-Shamji Industrial Estate, 9 LBS Marg, Kurla (W), Mumbai 400070
Tel: (022) 4215 5441 / 93204 22266

Hyderabad: H.NO-3-41-66, Sri Laxmi Nilayam, South Block, West Marredpally
Secunderabad 500026, Tel: 93468 32266

Lucknow: 4/11A, Ground Floor, Vibhav Khand, Gomti Nagar, Lucknow 226010
Tel: 93361 20941

E-mail: sales@nidekindia.com, Website: www.nidekindia.com